

**Instrument**

Vacuum manifold and separations system

**Model / Number**


Multiscreen

**Manufacturer**

Millipore

**Location**

Life Science 250H

**Instrument Description**

Vacuum manifold is designed for use with 96- and 384-well filter and receiver plates. This manifold allows rapid washing and/or collection of samples. The separations system is a filtration system designed to simplify all types of separations using a 96-well format, including biochemical assays, nucleic acid purifications, high throughput sample preparation, and sample preparation.

**User Charges**

Not applicable

**User Information**

User must have key card access to Life Science 250. User must be trained to use this equipment. The user is responsible for cleaning and performing any per use maintenance, and providing any disposable supplies required for operation. Any damage caused by the user is the responsibility of the user to repair or replace.