

Instrument

Centrifuge (x2)

Model / Number

Allegra X15-R
(Rotor SX4750)

Manufacturer

Beckman-Coulter

Location

Life Science 250 and Life
Science 250E

**Instrument Description**

Centrifuge with swinging bucket rotor. Specifications: Set speed 200 to 10,200 rpm in 10rpm increments; speed display-rpm or rcf: set time up to 99 hours 59 minutes or continuous hold; set temp -10 to 40°C in 1 degree increments; SX4750 rotor-max rpm 4,750; max rcf 5,250xg. Two buckets for 50mL tubes; two buckets for 15mL tubes.

User Charges

Not applicable

User Information

User must have key card access to Life Science 250. User must be trained to use this equipment. The user is responsible for balancing the load, cleaning and performing per run maintenance, completing the centrifuge log and providing any disposable supplies required for operation. Any damage caused by the user is the responsibility of the user to repair or replace.