

Instrument

Centrifuge

Model / Number

Avanti J-26XP (Rotor JA 25.5)

Manufacturer

Beckman-Coulter

Location

Life Science 250E

**Instrument Description**

Centrifuge max speed up to 82,000xg (26,000rpm). Temperature range -10 to 40°C. Rotor JA 25.5 max speed 25,000rpm (75,600xg); 8 x 50mL tubes (not conical tubes – must be Beckman tubes).

User Charges

Not applicable

User Information

User must have key card access to Life Science 250. User must be trained to use this equipment. The user is responsible for balancing the load, cleaning and performing per run maintenance, completing the centrifuge log and providing any disposable supplies required for operation. Any damage caused by the user is the responsibility of the user to repair or replace.