

Instrument
Centrifuge

Model / Number
5415R

Manufacturer
Eppendorf

Location
Life Science 250



Instrument Description

Up to 13,200 rpm / 16,100 xg. 24 x 1.5mL tubes. Fixed angle rotors.

User Charges

Not applicable

User Information

User must have key card access to Life Science 250. User must be trained to use this equipment. The user is responsible for balancing the load, cleaning and performing per run maintenance, completing the centrifuge log and providing any disposable supplies required for operation. Any damage caused by the user is the responsibility of the user to repair or replace.